

**Manual wheelchair
with electric verticalizer
User's manual**

WR6129F1



Contents

1. General description	1
2. Product performace	1
3. Product parameters	1
4. Structure and composition	1
5. Scope of application	2
6. Contraindications, precautions, warnings, and prompt instructions	2
6.1. Contraindications	2
6.2. Precautions, warnings and prompt instructions	2
6.3. Interpretation of graphics, symbols, abbreviations, etc. used in medical device labels	3
7. Installation	5
7.1. Installation	5
8. maintenance	7
9. Validity period	7
10. Method of cleaning, disinfecting, packaging and sterilizing	7
10.1 Basic guide	7
10.2 Manual cleaning and disinfection	7
11 Special storage, transportation conditions, methods	7
11.1. Storage	7
11.2. transport	7
12 Parts List	7
13 After sales and repair services	8

1. General description

This product is the Manual wheelchair with electric verticalizer, which produced according to the technical requirements of medical device products.

2. Product performance

Material: Powder coating steel chair frame

Flip-up armrest

Comfortable safety largest belt

Electric control stand up

Manual control the movement.

3. Product parameters(Unit mm)

Total width: 65cm

Seat width: 46cm

Rear wheel diameter: 60cm


Front wheel diameter: 20cm

Seat height: 58cm

Total height: 104-138cm

4 Structure and composition

This product is mainly composed of frame assembly, armrest assembly, front and rear wheel assembly, brake assembly, seat cushion, foot pedal etc.



5.Scope of application


This product is an important tool for rehabilitation. It is not only a means of transportation for physically handicapped and people with reduced mobility, but more importantly, it enables them to exercise and participate in social activities with the help of wheelchairs.

6 Contraindications, precautions, warnings, and prompt instructions

6.1 Contraindications: None

6.2 Precautions


- 1、 注意 The wheelchair has a rated load of 100Kg. When the user weighs more than 100Kg, please use it carefully or contact the manufacturer.



- 2、 注意 When the user pushes the wheelchair by himself, please keep driving at a constant speed and keep the speed of 3-5km/h..




- 3、 注意 When the wheelchair is driving on a sloping ground, it should be carried out under the condition that the user's physical condition is good, and the ground tilt angle is less than 10°. When the tilt angle is greater than 10°, it should be driven by others behind the back.



- 4、 警告 Users should not stand on the footrest during use to prevent the wheelchair from falling over.



5、 警告 When the wheelchair is in exercise, it is suitable to be carried out on a flat ground. When the current face encounters obstacles, it should be bypassed to avoid the risk of wheelchair falling.




6、 警告 When using a wheelchair on an inclined road, do not dump the wheelchair and suddenly change direction. When you are downhill, you cannot suddenly brake to prevent the risk of tipping over.



7、 警告 Before using the standing function, it is necessary to tie the safety straps of the upper and lower limbs, And make sure the handbrake is in the braking position

8、 警告 The hand controller is used to control the standing function of the wheelchair. This “” button controls standing up and this “” button controls standing down.

8. The light on the control box will not light up on its own. If want to check the battery level, you need to press the button "ON". If it shows green, it is fully charged, if it shows orange, it is about half charged, and if it shows red, it is almost no power.



The interpretation of the graphics, symbols, abbreviations used in the label, such as:

6.3 Interpretation of graphics, symbols, abbreviations, etc. used in medical device labels

symbols	Description
	Indicates a dangerous state of death or serious injury if not handled in accordance with the warnings.
	Indicates a dangerous state that could result in death or serious injury if not handled in accordance with the warnings. (1) Confirm that the specified power source is used. If a non-designated power source is used, fire or electric shock may occur.
	Indicates a dangerous state that may result in minor or moderate injury if not handled in accordance with the warning.

	Reference instructions
	note. Reference accompanying documents
	Avoid sun exposure
	Prevent heat and radiation sources
	Avoid rain
	Lower temperature limit
	Upper temperature limit
	Temperature limit
	Period of use
	Date of manufacture
	Batch code
	Classification number
	Serial number
	Temperature limit
note	Indicates a condition that may result in property damage if not operated as a warning.

7 Installation


7.1. Installation

Correct installation method


7.1.1 Open the package, take out the wheelchair backrest and wheelchair frame.




7.1.2 Loosen the backrest tube screws first.




7.1.3 Insert the backrest and tighten the screws.



7.1.4 Loosen the screw first, then insert the handle and tighten the screws.



7.1.4 This is the changer, when using, need to connect the Power cord and control box.



8 maintenance

Be sure to check the following items before use, so that the wheelchair is always in a safe state.:

8.1  注意 Check the screws on the universal wheel, drive wheel, brake and other parts. If it is loose, please lock it.

8.2  注意 Check the battery power regularly to ensure that the standing function can be used.

9 Validity period

1 year

10 Method of cleaning, disinfecting, packaging and sterilizing

10.1 Basic guide

Please comply with all relevant regulations on cleaning, disinfection and sterilization, as well as the internal system of the unit.

10.2 Manual cleaning and disinfection



Attention:

Wipe only with a damp cloth (wipe disinfection). It is forbidden to immerse the wheelchair in liquid.

Wipe the outer surface with a soft cloth dampened with water, mild soapy water or isopropyl alcohol

Can remove stubborn stains with a neutral detergent containing ammonia.

11 Do not use abrasive cleaners or solvents, as this may damage the surface or label

12 Special storage, transportation conditions, methods

11.1. Storage

Stored at a temperature of $-25^{\circ}\text{C} \sim 70^{\circ}\text{C}$, relative humidity 10% RH \sim 95% RH, atmospheric pressure: 700hPa \sim 1060hPa, no corrosive gas and well-ventilated clean room.

11.2. transportation

Products should be protected from heavy pressure, direct sunlight and rain when transported, and should be in accordance with the order contract.

Parts List

SN	Name	Number	remarks
1	wheelchair	1pcs	
2	Activity footrest	1pair	
3	Safety belt	1pcs	
4	Repair tool	1pcs	
5	Instruction Manual	1pcs	

12: After sales and repair services

From the date of your purchase of the company's wheelchair products; in normal use, if the wheelchair fails, you can get one year free warranty service. The company is responsible for free warranty or replacement. If there is a problem after one year, the company can be responsible for the repair, but reasonable recovery of the cost; the main frame of the company's products is guaranteed for 1 year.

For consumable parts, such as the normal wear and tear of large and small tires, it is not covered by the warranty.